















Piliers de cicatrisation

Marquage au laser pour une identification aisée du diamètre et de la hauteur



	Ø 3.4 mm			
				
Dimensions	D: Ø 3.4 mm H: 2 mm	D: Ø 3.4 mm H: 3 mm	D: Ø 3.4 mm H: 5 mm	D: Ø 3.4 mm H: 7 mm
Code	HSD3.4-2-CHC	HSD3.4-3-CHC	HSD3.4-5-CHC	HSD3.4-7-CHC
N° de Réf.	7311	7312	7313	7314

	Ø 3.8 mm			
				
Dimensions	D: Ø 3.8 mm H: 2 mm	D: Ø 3.8 mm H: 3 mm	D: Ø 3.8 mm H: 5 mm	D: Ø 3.8 mm H: 7 mm
Code	HSD3.8-2-CHC	HSD3.8-3-CHC	HSD3.8-5-CHC	HSD3.8-7-CHC
N° de Réf.	7315	7316	7317	7318

	Ø 4.2 mm			
				
Dimensions	D: Ø 4.2 mm H: 2 mm	D: Ø 4.2 mm H: 3 mm	D: Ø 4.2 mm H: 5 mm	D: Ø 4.2 mm H: 7 mm
Code	HSD4.2-2-CHC	HSD4.2-3-CHC	HSD4.2-5-CHC	HSD4.2-7-CHC
N° de Réf.	7319	7320	7321	7322